


TarryBoards™

TarryBoards are manufactured in our Connecticut facility for hospitals all over the world. Sizes range from Micro-Preemie (¾" x 2") to Adult X-Large (4" x 17").


- Provides comfort for the patient while supporting the I.V. site
- Handy adjustment tab at each end
- Lined with a thick foam pad to conform gently to the limb
- Covered in soft spun fabric
- Single patient use
- Latex-free
- Made in the U.S.A.



Preemie and Micro-Preemie Armboards



Preemie, Neonatal, Pediatric and Child-Size Armboards



Adult Armboards

Description		Item#	Qty/Box	Sizes
Micro-Preemie	Semi-Flexible	T-52839	50/100	¾" x 2"
	Unbreakable	T-92639	50/100	
Preemie	Semi-flexible	T-52840	50/100	1" x 2½"
	Unbreakable	T-92640	50/100	
	Ultra-unbreakable	T-92670	50/100	
Preemie Plus	Semi-flexible	T-52831	50/100	1½" x 3"
	Unbreakable	T-92631	50/100	
	Ultra-unbreakable	T-92671	50/100	
Neonatal	Semi-flexible	T-52841	50/100	1½" x 4"
	Unbreakable	T-92641	50/100	
	Ultra-unbreakable	T-92677	50/100	
Neonatal Plus	Semi-flexible	T-52832	50/100	2" x 4"
	Unbreakable	T-92632	50/100	
	Ultra-unbreakable	T-92672	50/100	
Pediatric	Semi-flexible	T-52842	50/100	2" x 5"
	Unbreakable	T-92642	50/100	
	Ultra-unbreakable	T-92676	50/100	
Pediatric Plus	Unbreakable	T-92662	50/100	2" x 6"
	Ultra-unbreakable	T-92673	50/100	
Pediatric Plus Wide	Unbreakable	T-92664	50/100	2½" x 6"
	Ultra-unbreakable	T-92674	50/100	
Child Narrow	Ultra-Unbreakable	T-1585	50/100	1½" x 8½"
Child	Unbreakable	T-52845	50/100	3" x 7"
	Ultra-unbreakable	T-92678	50/100	
Child Plus	Unbreakable	T-52893	50/100	3" x 9"
	Ultra-unbreakable	T-92675	50/100	
Adult		T-52843	50/100	4" x 9"
Adult Large		T-52844	50/100	4" x 14"
Adult X-Large		T-52846	50/100	4" x 17"

Manufactured exclusively by
 Tarry Manufacturing
 Samples available upon request.

customerservice@tarrymfg.com

www.tarrymfg.com