

TarryGels™

TarryGels are soft gel filled pillows designed to provide a comfortable cushion to meet the support needs of micro-preemie, preemie, and full-term infants. These pillows can help reduce the need for frequent re-positioning or to make immobilized babies more comfortable. Each TarryGel comes with a fully washable, attractive, flame retardant cover. The covers are also available for purchase separately.

USES:

- Head pillow
- Mattress for micro-preemies
- Positioning Aid or Prone Roll
- Single patient, hospital use only

PRODUCT DESCRIPTION:

- Pillows are made from soft, DEHP-free polyethylene film that does not react with the gel fill. The film also prevents the gel fill from evaporating
- Pillow covers are a flame retardant polyester blend material that can be washed and reused (single patient only)
- Gel is a non-toxic, water based food gel designed to be resistant to bacterial and fungal growth
- Available in 3 sizes: 3" x 6", 6" x 6", and 6" x 9"
- Latex-Free
- Made in the USA

PRECAUTIONS/DIRECTIONS:

- **Warning!** *Danger of suffocation!* Never leave baby face down and unattended.
- **Warning!** Always closely monitor the baby
- **Caution!** Do not x-ray as small gel air pockets may distort x-ray results
- **Caution!** Do not microwave as hot spots may form
- Warm before using
- Always use a pillow cover



Description	Item#	Qty/Box	Sizes
Small Gel Pillow,	T 92332	Box of 12	3" x 6"
Medium Gel Pillow	T 92335	Box of 12	6" x 6"
Large Gel Pillow	T 92338	Box of 12	6" x 9"
Small Gel Pillow Covers	T 92332 C	Box of 50	Fits 3" x 6" pillows
Medium Gel Pillow Covers	T 92335 C	Box of 50	Fits 6" x 6" pillows
Large Gel Pillow Covers	T 92338 C	Box of 50	Fits 6" x 9" pillows

Manufactured exclusively by
 Tarry Manufacturing
 Samples available upon request.

customerservice@tarrymfg.com

www.tarrymfg.com