

# BD Enteral Syringe


Designed to meet your needs.



Helping all people  
live healthy lives

**BD Medical**  
1 Becton Drive  
Franklin Lakes, NJ 07417

**BD Enteral Syringe — 1 mL – 60 mL**



**U.S. MANUFACTURED  
SOURCE**

**FDA 510(K)  
CLEARANCE**


**VISIBLE DIFFERENTIATION  
FROM PARENTERAL SYRINGE**



## We put everything you asked for into our Enteral Syringe.

We're helping you safely deliver nourishment and oral medications.

And there's no need to pay for costly syringe pump software upgrades. BD's Enteral Syringe does it all – and more.



### 1. Amber-colored plunger rod

eases oral delivery in low-light conditions and provides contrast between medicine and barrel

### 2. Clear barrel

provides simple view of medication in any situation

### 3. Bold print

makes reading scale markings easier

### 4. Non-luer compatible BD UniVia™ connection

designed to reduce risk of tubing set misconnection




### Other features

- **U.S. made Enteral Syringe**  
with FDA 510(K) clearance
- **Sterile and packaged individually**
- **Quality performance**  
equivalent to current BD Parenteral Syringes
- **Differentiation in design**  
based on clinical feedback

**No need for costly pump software upgrades.  
All at a cost you'll find hard to pass up.**

# The right solution for safety and accuracy.

Our Enteral Syringe will deliver on all expectations.



Amber-colored plunger rod

Non-luer fitting  
BD UniVia™ connection

BD Catalog No.	Product Description	Shelf Carton Quantity	Case Quantity
305851	1 mL Syringe	200 each	1,600 each or 8 shelf cartons
305853	3 mL Syringe	200 each	800 each or 4 shelf cartons
305855	5 mL Syringe	125 each	500 each or 4 shelf cartons
305857	10 mL Syringe	100 each	400 each or 4 shelf cartons
305859	20 mL Syringe	40 each	160 each or 4 shelf cartons
305861	30 mL Syringe	40 each	160 each or 4 shelf cartons
305863	60 mL Syringe	40 each	160 each or 4 shelf cartons
302435	Oral Tip Cap	100 each	1,000 each or 10 shelf cartons



**The BD Enteral Syringe provides syringe pump performance,** sterility and differentiation of design and appearance – in one solution.

Think BD for Enteral delivery. BD has mastered delivering products that perform safely and accurately.

**To learn more or to request samples, contact your local BD Medical Sales Consultant.**

**Or contact us at (201) 847-5411**  
**[Hypodermic\\_Tech\\_Support@bd.com](mailto:Hypodermic_Tech_Support@bd.com)**

**Your BD blueprint.**  
Constructing something better.



**BD Medical**  
1 Becton Drive  
Franklin Lakes, NJ 07417  
[www.bd.com](http://www.bd.com)