




Peripherally Inserted Central Catheter with Peel-Away Introducer Sheaths



The hub is clearly marked with "PICC" and French size to readily identify catheter type.

The Footprint® introducer sheaths have unique design features to enhance first stick success.



Peripherally Inserted Central Catheter

- Designed especially for neonatal/pediatric patients
- Silicone or polyurethane catheter with integrated extension set
- Soft stabilization wing for enhanced securement
- Indelible markings every centimeter
- Fully radiopaque
- Latex and DEHP free

Introducer Sheath

- Notched needle with minimal lie distance for immediate visualization of flashback designed to increase first stick success
- Low profile, compact introducer designed to enhance small vein access



Footprint[®]
medical inc.

210.226.2600

www.footprintmed.com



Ordering Information for:

Peripherally Inserted Central Catheters

Silicone Peripherally Inserted Central Catheters

S1PIC1.9-S	Silicone PICC Complete Kit <i>(with catheter & introducer)</i>	1.9 Fr.	5/cs
S1PIC1.9-C	Silicone PICC Catheter Only (28cm)	1.9 Fr.	5/cs
S1PIC1.9-L	Silicone PICC Catheter Only (43cm)	1.9 Fr.	5/cs

Polyurethane Peripherally Inserted Central Catheters

P1PIC1.4-S	Polyurethane PICC Complete Kit <i>(with catheter & introducer)</i>	1.4 Fr.	5/cs
P1PIC1.4-C	Polyurethane PICC Catheter Only (28cm)	1.4 Fr.	5/cs
P1PIC1.9-S	Polyurethane PICC Complete Kit <i>(with catheter & introducer)</i>	1.9 Fr.	5/cs
P1PIC1.9-C	Polyurethane PICC Catheter Only (28cm)	1.9 Fr.	5/cs

Peel-Away Introducer Sheaths

SHEATH1.4	PICC Sheath Introducer Only	1.4 Fr.	5/cs
SHEATH1.9	PICC Sheath Introducer Only	1.9 Fr.	5/cs

"Helping tiny babies make their own footprints in the world..."

Distributed By:



Phone: 847.891.6210
Fax: 847.891.6215
www.accessmedicalsupply.com

Access Medical Supply, Inc.



12727 Cimarron Path
San Antonio, TX 78249
210.226.2600
www.footprintmed.com