



Silicone Enteral Feeding Tube



The open end design, with precision drilled side holes, enhances fluid flow and decreases potential feeding stagnation and mucosal irritation.

The hub is clearly marked with "ORAL" and French size to readily identify the tube as an enteral feeding tube.



Silicone Enteral Feeding Tube

- Made of soft, biocompatible medical grade silicone designed to enhance extended dwell protocols
- Labeling supports ISMP recommendations that all distal ports and tubing of access lines clearly indicate what is the use of the port/line
- Bold, indelible markings every centimeter
- Color coded for patient safety
- Enteral connector designed to prevent misconnections
- Integrated closure plug for easy cleaning
- Radiopaque stripe for confirmation of tip placement by x-ray
- Latex and DEHP free



Ordering Information for:

Silicone Feeding Tubes

with Enteral Connector



Silicone 40cm Enteral Feeding Tubes

S1FTS5.0-EO	Silicone Feeding Tubes	5.0 Fr.	10/cs
S1FTS6.5-EO	Silicone Feeding Tubes	6.5 Fr.	10/cs
S1FTS8.0-EO	Silicone Feeding Tubes	8.0 Fr.	10/cs

Silicone 60cm Enteral Feeding Tubes

S1FTM5.0-EO	Silicone Feeding Tubes	5.0 Fr.	10/cs
S1FTM6.5-EO	Silicone Feeding Tubes	6.5 Fr.	10/cs
S1FTM8.0-EO	Silicone Feeding Tubes	8.0 Fr.	10/cs

Silicone 90cm Enteral Feeding Tubes

S1FTL5.0-EO	Silicone Feeding Tubes	5.0 Fr.	10/cs
S1FTL6.5-EO	Silicone Feeding Tubes	6.5 Fr.	10/cs
S1FTL8.0-EO	Silicone Feeding Tubes	8.0 Fr.	10/cs

"Helping tiny babies make their own footprints in the world..."

Distributed By:



Phone: 847.891.6210
Fax: 847.891.6215
www.accessmedicalsupply.com

Access Medical Supply, Inc.



12727 Cimarron Path
San Antonio, TX 78249
210.226.2600
www.footprintmed.com