



Urinary Drainage Catheter



The 100ml burette provides accurate measurement of urinary output.

The hub is clearly marked with "URINE" and French size to readily identify the catheter type.



Urinary Drainage Catheter

- Made of soft, biocompatible medical grade silicone designed to enhance extended dwell protocols
- Bold, indelible markings every centimeter
- Catheter has a rounded tip with two precision drilled side holes to reduce trauma during insertion and minimize urethral irritation
- Luer-Lok™ hub with integrated closure plug
- Latex and DEHP free with radiopaque stripe
- Also available in Urine Collection Kit
- Closed System Kit Contains:
 - Silicone urinary catheter
 - 100ml Burette
 - Antiseptic swab sticks
 - Sterile catheter lubricant
 - Fenestrated drape
 - Poly-lined drape
 - Gauze pads



Ordering Information for:

Urinary Drainage Catheter



Silicone Urinary Drainage Catheters

S1UDC3.5	3.5 Silicone Urinary Drainage Catheter	10/cs
S1UDC5.0	5.0 Silicone Urinary Drainage Catheter	10/cs
S1UDC6.5	6.5 Silicone Urinary Drainage Catheter	10/cs
S1UDC8.0	8.0 Silicone Urinary Drainage Catheter	10/cs

Silicone Urinary Drainage Kit without Catheter

UDKNC	Closed System Urinary Collection Kit without Catheter	5/cs
-------	---	------

Silicone Urinary Collection Kit with Catheter

S1UDK3.5	3.5 Fr. Silicone Urinary Collection Kit with Catheter	5/cs
S1UDK5.0	5.0 Fr. Silicone Urinary Collection Kit with Catheter	5/cs
S1UDK6.5	6.5 Fr. Silicone Urinary Collection Kit with Catheter	5/cs
S1UDK8.0	8.0 Fr. Silicone Urinary Collection Kit with Catheter	5/cs

"Helping tiny babies make their own footprints in the world..."

Distributed By:



Access Medical Supply, Inc.

Phone: 847.891.6210
Fax: 847.891.6215
www.accessmedicalsupply.com



12727 Cimarron Path
San Antonio, TX 78249
210.226.2600
www.footprintmed.com