

EyeMax2™

- New easy-to-use design
- Adjustable band for a perfect fit
- Bigger eye pad to reduce pressure and increase UV protection



Phototherapy Masks

Maxtec® is pleased to introduce the all new EyeMax™2 phototherapy mask. Based on valuable feedback from our customers, we have improved our design to assure optimum comfort and protection for our patients during phototherapy treatment. The EyeMax™2 is designed to adapt to the baby's changing head shape. The single-piece, wrap-around design attaches at two points on either side of the baby's head and can be independently adjusted to prevent unwanted movement and allow a perfect fit. Our special UV light eye protection pad has been ergonomically designed to prevent any "light leaks". The EyeMax™2 is available in three sizes and come in packs of 20.

Specifications

Regular	12.6" - 14.9" (Occipital Frontal Circumference)
Preemie	10.4" - 12.6" (Occipital Frontal Circumference)
Micro	7.87" - 10.4" (Occipital Frontal Circumference)

Part Numbers

Regular (comes in a pack of 20)	R300P01	Preemie (comes in a pack of 20)	R300P02
Micro (comes in a pack of 20)	R300P03		