



TORTLE NICU MIDLINER

Precision Head Positioning System for micro-preemies & preemies.

MidLiner Benefits:

Supports Supine Midline Position

- The AAP recommends supine midline position for micro-preemies during the first 72+ hours of life to help prevent intra-ventricular hemorrhage (IVH)

Velcro front opening and adjustable tabs

- Simple and easy to apply which minimizes handling and stress on baby
- Provides easy access for scalp IVs and head ultrasounds
- Ensures a proper and comfortable fit
- Reduces staff intervention time

Cannula tabs and numerous device attachments

- Securely holds most nasal CPAP apparatuses, nasal cannulas including RAM, as well as feeding tubes
- Includes additional straps to hold most bilishades
- X-ray friendly and compatible with ventilator

Simplifies transition to all positions, including supine midline, side lying, and prone (log rolling)

- No need to detach and re-apply respiratory or feeding apparatus
- Helps prevent and treat Dolichocephaly/Scaphocephaly

Available in three sizes:

SIZE	CASE#	WEIGHT	HEAD CIRCUMFERENCE
ELBW	65911	.50-1.25 Kg	24-27 cm
MICRO	65901	1.25-2.00 Kg	27-32 cm
PREMIE	65903	2.00-3.50 Kg	32-37 cm

10 units per case. Machine washable.




Supports AAP Back-to-Sleep Campaign



The FDA-cleared Turtle NICU MidLiner is a simple and efficient head repositioning system that assists staff in providing optimal patient care.

Three Primary Positions

1 - Supine Midline



2 - Side Lying



3 - Prone



Distributed By:



Access Medical Supply, Inc.

Phone: 847.891.6210
Fax: 847.891.6215
www.accessmedicalsupply.com