

SKINTACT Cool Contact • Electrosurgical Grounding Plates

"SKINTACT Cool Contact electro-surgical grounding plates provide safety of the highest degree. They are compatible with all makes of generators on the market and are totally compliant with AAMI electro-surgical standards guaranteeing patient safety and your peace of mind."


















SKINTACT Cool Contact • Electrosurgical Grounding Plates

The professional's choice

- ❑ Unique C-LINE manufacturing produces ultra soft hydrogel which conforms perfectly to the body contours providing maximum safety.
- ❑ Low skin contact impedance ensures the plate stays as cool as possible during use.
- ❑ Flexible low profile foam, easy to mould provides excellent security.
- ❑ An enhanced adhesive border provides extra security should the plate be pulled inadvertently or come into contact with prepping solutions.
- ❑ Available in a comprehensive range of sizes, adult/pediatric, single/split and with or without pre-attached cable.

SKINTACT Cool Contact Dispersive Grounding Plates

					
Order Code	RS05	RS06	RS25	RS26	RS27
Type	single	single	split	split	split
Size	adult	pediatric	adult	pediatric	adult
Conductive Area	141 cm ²	76 cm ²	125 cm ²	68 cm ²	131 cm ²
Packaging Standard	5 / pouch 50 / box	1 / pouch 50 / box	5 / pouch 50 / box	1 / pouch 50 / box	5 / pouch 50 / box
with 3 m pre-attached disposable cable					
Order Code	RS65B	RS66B	RS35A	RS36A	RS37A
Connector	Non-REM	Non-REM	REM	REM	REM
Packaging Standard	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box
with 3 m pre-attached disposable cable					
Order Code	RS65D	RS66D	RS35D	RS36D	RS37D
Connector	6.3 mm jack	6.3 mm jack	6.3 mm jack	6.3 mm jack	6.3 mm jack
Packaging Standard	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box	1 / pouch 50 / box

SKINTACT Cool Contact Reusable Dispersive Cables and Adaptors

							
Order Code	NEK300A	NEK300B	NEK300D	NEK300E	NEK300F	NEK300G	NEK300H
Connector (compatible)	REM (Valleylab)	Non-REM (Valleylab)	6.3 mm jack (standard)	4 mm banana	8 Pin (Siemens)	2 Pin adjustable (Eschmann)	6.3 mm jack (Erbe)
Cable	3 m	3 m	3 m	3 m	3 m	3 m	3 m
Order Code	NEK500A	NEK500B	NEK500D	NEK500E	NEK500F	NEK500G	NEK500H
Cable	5 m	5 m	5 m	5 m	5 m	5 m	5 m
Order Code		NEKA-AD		NEKA-BE		RAZOR	
Description		Adaptor A/B-D		Adaptor B-E		Prep razor	

SKINTACT Cool Contact Dispersive Grounding Plates, Cables and Adaptors meet the ANSI/AAMI standard HF18:2001 for Electrosurgical Devices.