

## Table of Contents

### **NeoDevices**

Umbilical Catheterization Trays	Page 1
Forceps	Page 1
PICC Insertion Tray	Page 1
Feeding Tubes (Silicone)	Page 2
Feeding Tubes (Polyurethane)	Page 2
Urinary Drainage Catheters (Silicone)	Page 2
Closed Urinary Collection Kits (Silicone)	Page 3
Neonatal Lumbar Puncture Needles	Page 3
Neonatal Lumbar Puncture Kits	Page 3
Prewired Hydrogel Electrodes	Page 4
Skin Temperature Probe	Page 4

### **Infusion Devices**

Enteral Feeding Extension Sets - Orange Markings	Page 5
Syringe Pump Tubing Extension Sets	Page 5-11

**NeoDevices**

**NeoDevices Umbilical Catheterization Tray**



Description	Size	Reorder No.	Units/Case	Cross Ref.	Cross Ref.
Individually packaged forceps	N/A	CT	5/case	CATHTRAY	4070007

**NeoDevices Forceps**



Description	Size	Reorder No.	Units/Case	Cross Ref.
Forceps	N/A	F275	10/case	NEOCEPS

**NeoDevices PICC Insertion Tray**



Description	Size	Reorder No.	Units/Case
PICC Insertion Tray without catheter	N/A	PIT	5/case

## NeoDevices

### NeoDevices Feeding Tubes - Silicone



Description	Size	Reorder No.	Units/Case	Cross Ref.	Cross Ref.
Feeding Tube – short (40cm length)	5.0 fr	FTS5.0S	10/case	S1FTS5.0	4165019
Feeding Tube – short (40cm length)	6.5fr	FTS6.5S	10/case	S1FTS6.5	4166519
Feeding Tube – short (40cm length)	8.0 fr	FTS8.0S	10/case	S1FTS8.0	4168019
Feeding Tube – medium (60cm length)	5.0 fr	FTM5.0S	10/case	S1FTM5.0	4165029
Feeding Tube – medium (60cm length)	6.5fr	FTM6.5S	10/case	S1FTM6.5	4166529
Feeding Tube – medium (60cm length)	8.0 fr	FTM8.0S	10/case	S1FTM8.0	4168029
Feeding Tube – long (90cm length)	5.0 fr	FTL5.0S	10/case	S1FTL5.0	4165039
Feeding Tube – long (90cm length)	6.5 fr	FTL6.5S	10/case	S1FTL6.5	4166539
Feeding Tube - long (90cm length)	8.0 fr	FTL8.0S	10/case	S1FTL8.0	4168039

### NeoDevices Polyurethane feeding tubes



Description	Size	Reorder No.	Units/Case
Feeding Tube – short (40cm length)	4.0 fr	FTS4.0P	10/case
Feeding Tube – short (40cm length)	5.0 fr	FTS5.0P	10/case
Feeding Tube – short (40cm length)	6.5 fr	FTS6.5P	10/case
Feeding Tube – short (40cm length)	8.0 fr	FTS8.0P	10/case
Feeding Tube – medium (60cm length)	4.0 fr	FTM4.0P	10/case
Feeding Tube – medium (60cm length)	5.0 fr	FTM5.0P	10/case
Feeding Tube – medium (60cm length)	6.5 fr	FTM6.5P	10/case
Feeding Tube – medium (60cm length)	8.0 fr	FTM8.0P	10/case
Feeding Tube – long (90cm length)	5.0fr	FTL5.0P	10/case
Feeding Tube – long (90cm length)	6.5 fr	FTL6.5P	10/case
Feeding Tube - long (90cm length)	8.0 fr	FTL8.0P	10/case

### NeoDevices Urinary Drainage Catheter - Silicone



Description	Size	Reorder No.	Units/Case	Cross Ref.	Cross Ref.
Urinary Drainage Catheter	3.5 fr	UC3.5	10/case	S1UDC3.5	None
Urinary Drainage Catheter	5.0 fr	UC5.0	10/case	S1UDC5.0	4195005
Urinary Drainage Catheter	6.5 fr	UC6.5	10/case	S1UDC6.5	None
Urinary Drainage Catheter	8.0 fr	UC8.0	10/case	S1UDC8.0	4198005


## NeoDevices

### NeoDevices Closed Urinary Collection Kits - Silicone



Description	Size	Reorder No.	Units/Case	Cross Ref.	Cross Ref.
Urinary Collection Kit	3.5 fr	UK3.5	5/case	S1UDK3.5	None
Urinary Collection Kit	5.0 fr	UK5.0	5/case	S1UDK5.0	4195007
Urinary Collection Kit	6.5 fr	UK6.5	5/case	S1UDK6.5	None
Urinary Collection Kit	8.0 fr	UK8.0	5/case	S1UDK8.0	4198007

### NeoDevices Neonatal Lumbar Puncture Needles



Description	Size	Reorder No.	Units/Case	Cross Ref.	Cross Ref.
22 Gauge 1 inch Spinal Needle	22 ga	LN2210	10/case	LPN-22-1.0	4011019
22 Gauge 1½ inch Spinal Needle	22 ga	LN2215	10/case	LPN-22-1.5	4011529


### NeoCare Neonatal Lumbar Puncture Kits



Description	Size	Reorder No.	Units/Case	Cross Ref.
Kit w/22 Gauge 1 inch Spinal Needle	22 ga	LPK-22-1.0	5/case	4011019
Kit w/25 Gauge 1 inch Spinal Needle	25 ga	LPK-25-1.0	5/case	None

## NeoDevices

### Prewired Hydrogel Electrodes



Description	Size	Reorder No	Units/Case
Premie Prewired Hydrogel Electrodes	7/8 X 7/8	E-PC1-6C	100 x 3/pouch
Newborn Prewired Hydrogel Electrodes	7/8 X 1 1/4	E-PC2-6C	100 x 3/pouch
Pediatric Prewired Hydrogel Electrodes	1 1/4 X 1 1/4	E-PC3-6C	100 x 3/pouch

### Skin Temperature Probe Covers



Description	Reorder No.	Units Per
"Cuddle Covers" Hydrogel Probe Cover	CC-100	100/box
Heart Shaped Hydrogel Cover	HH-100	100/box
"Cuddle Covers" Foam Probe Cover	FC-100	100/box

## INFUSION DEVICES

### Enteral Feeding Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)

IS60EN



IS60ENY-C



IS60ENSC



IS60ENY




Description	Size	Reorder No.	Units/Case
Enteral tubing, <b>Orange Stripe</b> (2.0ml priming volume) Female luer lock and <b>Orange</b> male <b>enteral</b> adapter (non-IV)	60 inch	IS60EN	50/case
Enteral tubing, <b>Orange Stripe</b> w/ Y ext. site(2.0ml priming volume) Female luer lock and Orange male enteral adapter (non-IV)	60 inch	IS60ENY-C	50/case
Enteral tubing, <b>Orange Stripe</b> w/ stopcock(2.0ml priming volume) Female luer lock and Orange male enteral adapter (non-IV)	60 inch	IS60ENSC	50/case
Enteral tubing, <b>Orange Stripe</b> w/ Y site (2.0ml priming volume) Female luer lock and Orange male enteral adapter (non-IV)	60 inch	IS60ENY	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
Childrens Med Venture	95017-A	60" Enteral Feeding Set	2.0 ml

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Microbore Extension Set (0.3ml priming volume) (microbore tubing for anesthesia/syringe/ambulatory pumps with female/male luer locks, & slide clamp)	60 inch	IS6003	50/case


#### Cross Reference

Company	Catalog #	Description	Priming Volume
Baxter	1M8511	60" Microbore Extension Set	0.40ml
Burron	ET-60M	60" Microbore Extension Set	0.30ml
Churchill	AMS-467	60" Microbore Extension Set	0.40ml
Medex	53-60-40	60" Microbore Extension Set	0.30ml
Medex	53-60-40C	60" Microbore Extension Set	0.30ml
Medlon	BC-576	60" Microbore Extension Set	0.30ml
MPS	EA06001	60" Microbore Extension Set	0.30ml

## INFUSION DEVICES

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Minibore Extension Set (1.0ml priming volume) (minibore tubing for anesthesia/syringe/ambulatory pumps with female/male luer locks & slide clamp)	60 inch	IS6010	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
Medex	53-60-35	60" Minibore Extension Set	1.0ml
MPS	EG06001	60" Minibore Extension Set	1.0ml

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Minibore Extension Set w/Stopcock (1.0ml priming volume) (minibore tubing for anesthesia/syringe/ambulatory pumps with female/male luer locks & slide clamp)	60 inch	IS6010SC	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
Medex		60" Minibore Extension Set	1.0ml

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Minibore Extension Set (1.8ml priming volume) (minibore tubing for syringe/ambulatory pumps with female/male luer locks & pinch clamp)	60 inch	IS6018	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918M160	60" Minibore Extension Set	1.7ml
Baxter	1M8521	60" Minibore Extension Set	1.9ml
Burron	ET-112	60" Minibore Extension Set	2.0ml
Churchill		60" Minibore Extension Set	
Medex	53-60-225	60" Minibore Extension Set	2.5ml
Medlon	BC-566	60" Minibore Extension Set	1.7ml
Medlon	BC-568	60" Minibore Extension Set	2.4ml
MPS	EB06001	60" Minibore Extension Set	1.8ml

## INFUSION DEVICES

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Dark Tubing Minibore Extension Set (1.8ml priming volume) (dark tubing for light sensitive fluids with female/male luer locks & pinch clamp)	60 inch	IS6018DT	50/case

### Syringe Pump Tubing Extension Set


(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Microbore Extension Set (0.2ml priming volume) (microbore tubing for anesthesia/syringe/ambulator pumps with female/male luer locks & slide clamp)	7 inch	IS0702	50/case

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Minibore Extension Set (2.0ml priming volume) (minibore tubing with anti-siphon valve, removable slide clamps, & male luer lock)	72 inch	IS72V	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918M172	72" Microbore Extension Set	2.1ml
AMD	9918M172S	72" Microbore Extension Set	2.1ml
AMD	9918M172V	72" Microbore Extension Set	2.1ml
Baxter/Bard	6466070	72" Microbore Extension Set	1.5ml
Baxter/Bard	6466120	72" Microbore Extension Set	2.8ml
Baxter/Bard	6466090	72" Microbore Extension Set	1.6ml
MPS	EB07201	72" Microbore Extension Set	2.0ml
MPS	EB07202	72" Microbore Extension Set	2.0ml



## INFUSION DEVICES

### Syringe Pump Tubing Extension Set


(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Minibore Extension Set (2.0ml priming volume) (minibore tubing for anesthesia/syringe/ambulatory pumps with female/male luer locks & slide clamp)	72 inch	IS7220	50/case

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



(for viscous medications/solutions)

Description	Size	Reorder No.	Units/Case
Macrobores Extension Set (5.0ml priming volume) (Kink resistant tubing, female luer lock, slide clamp, male luer lock)	30 inch	IS3050	50/case

### Syringe Pump Tubing Extension Set


(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Microbores Extension Set (0.7ml priming volume) (microbores tubing with .22 micron filter)	60 inch	IS60F	50/case

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Microbores Extension Set (0.18ml priming volume) (microbores tubing for anesthesia/syringe/ambulator pumps with female/male luer locks & slide clamp)	36 inch	IS3618	50/case


#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918M036	36" Microbores Extension Set	0.18ml
Baxter	1M8510	36" Microbores Extension Set	0.30ml
Burron	ET-36M	36" Microbores Extension Set	0.18ml
Churchill		36" Microbores Extension Set	
Marquette	MM-36	36" Microbores Extension Set	0.20ml
Medlon	53-36-40	36" Microbores Extension Set	0.20ml
Medlon	BC-574	36" Microbores Extension Set	0.18ml
MPS	EA03601	36" Microbores Extension Set	0.18ml

## INFUSION DEVICES

### Syringe Pump Tubing Extension Set


(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Minibore Extension Set (1.5ml priming volume) (female luer lock, slide clamp, male luer lock)	96 inch	IS9615	50/case

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Bifurcated Extension Set (2.1ml priming volume) (TPN/LIPID/PCA Combination with kink-resistant tubing, pinch clamps, and female/male luer locks)	8 inch	ISB21	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918MS2W5	2 Way 5" Extension Set	
Burron	ET-06Y	2 Way 6" Extension Set	1.4ml
Medex	BF-100	2 Way 6" Extension Set	0.6ml
Medlon	BC-588	2 Way 5" Extension Set	0.6ml
MPS	BC-00603	2 Way 6" Extension Set	0.85ml
MPS	BC-00602	2 Way 6" Extension Set	0.85ml

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Bifurcated Extension Set (0.3ml priming volume) (minibore tubing for NICU/PICU with slide clamps, female luer locks and male rotating hub)	4 inch	ISM2	50/case


#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918MI2W3	2 Way 3" Extension Set	0.16ml
AMD	9918MI2W6	2 Way 6" Extension Set	0.24ml
Baxter	2N1191	2 Way 4" Extension Set	
Burron	ET-103	2 Way 7" Extension Set	0.40ml
Medex	MBF-27	2 Way 4.5" Extension Set	0.30ml
Medex	MBF-30	2 Way 6" Extension Set	0.40ml
Medlon	BC-205	2 Way 5" Extension Set	0.20ml
Burron	ET-103	2 Way 7" Extension Set	0.40ml
Medex	MBF-27	2 Way 4.5" Extension Set	0.30ml
Medex	MBF-30	2 Way 6" Extension Set	0.40ml
Medlon	BC-205	2 Way 5" Extension Set	0.20ml
Medex	MBF-30	2 Way 6" Extension Set	0.40ml
Medlon	BC-205	2 Way 5" Extension Set	0.20ml

## INFUSION DEVICES

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)




Description	Size	Reorder No.	Units/Case
Y Extension Set (1.5ml priming volume) and female/male luer locks	8 inch	ISSBC21	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918MI2W5V	2 Way 5' Extension Set w/1 Valve	
AMD	9918MI2W2V6	2 Way 6" Extension Set w/2 Valves	
Baxter	1C8209	2 Way 6" Extension Set w/2 Valves	
Medex	BBF-26	2Way 5.5" Extension Set w/1 valve	0.50ml
Medlon	BC-200	2Way 3" Extension Set w/2 valves	0.15ml
Medlon	BC-211	2Way 3" Extension Set w/2 valves	

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Trifurcated Set (2.5ml priming volume) (kink-resistant tubing for ICU/PICU with color coded pinch clamps and female/male luer locks)	6 inch	IST31	50/case


#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918MS3W5	3 Way 6' Extension Set	
Medex	TF-55	3 Way 7" Extension Set	0.85ml
Medex	BC-587	3 Way 8" Extension Set w/Inj. Site	2.50ml
MPS	BE-00601	3 Way 6" Extension Set	1.50ml
MPS	BE-00611	3 Way Extension Set w/Inj. Site	1.80ml

## INFUSION DEVICES

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)




Description	Size	Reorder No.	Units/Case
Trifurcated Set (0.6 priming volume) (Minibore tubing for NICU/PICU with female/male luer locks and slide clamps)	8 inch	ISMT31	50/case

#### Cross Reference

Company	Catalog #	Description	Priming Volume
AMD	9918MI3W9	3 Way 9" Extension Set	0.40ml
AMD	9918MI3W7	3 Way 7" Extension Set	0.25ml
AMD	9918MI3W5	3 Way 5" Extension Set	0.25ml
Medex	MTF-37	3 Way 4" Extension Set	0.40ml
Medlon	BC-225	3 Way 7" Extension Set	0.50ml
Medlon	BC-559	3 Way 8" Extension Set	0.50ml
MPS	BG-00601	3 Way 6" Extension Set	0.40ml

### Syringe Pump Tubing Extension Set

(Non-DEHP PVC, No Latex, Lipid resistant)



Description	Size	Reorder No.	Units/Case
Mini Trifurcated Set (.33ml Priming Volume) (minibore tubing, open double T-connector, male luer lock, female luer locks, removable slide clamps)	2 inch	B1025	50/case